

MAKE SENSE OF KERATOCONUS

FIT SOFT | FIT HEALTHY | FIT THIN

KERA  SOFT THIN

A REVOLUTIONARY ULTRA-THIN SOFT LENS DESIGN FOR PEOPLE WITH KERATOCONUS AND IRREGULAR CORNEAS.

KeraSoft® Thin soft contact lenses feature a sophisticated optical design that delivers optimal visual acuity to the irregular cornea patient. This patented design minimizes overall thickness to maximize oxygen transmission and comfort while conforming to the shape of any cornea.

At a glance:

- A superb alternative to rigid corneal & scleral lenses
- Designed to fit a wide range of irregular corneas including all stages of keratoconus
- Comes in both spherical and toric parameters
- Highly customizable from every design aspect
- Ultra-thin design delivers enhanced peripheral fit

Custom Contact Lenses
METRO
OPTICS
800-223-1858
metro-optics.com

Silicone hydrogel for health, comfort, and excellent visual acuity for difficult-to-fit corneas.



Ease of fit

- A small, eight-lens diagnostic fitting set that aids in determining the correct lens parameters
- Pragmatic fitting procedure that quickly assesses lens fitting options
- Proven troubleshooting procedures for enhanced fitting accuracy

Patented design

- Front surface asphere or aspheric toric prism ballasted with balanced overall thickness and spherical aberration control
- When required, peripheries can be steepened or flattened independently of base curve
- Up to two sectors of the periphery can be modified independently with our proprietary Sector Management Control™

Comfort

- Silicone hydrogel — longer wearing comfort
- Provides high level of oxygen for crisp, clear, comfortable vision
- More oxygen compared to traditional hydrogel lenses

Standard Fitting Set Parameters

Base Curve	Diameter	Periphery	Power
7.80mm	14.5mm	STD	Plano
8.00mm	14.5mm	STD	Plano
8.20mm	14.5mm	STD	Plano
8.40mm	14.5mm	STD	Plano
8.60mm	14.5mm	STD	Plano
8.80mm	14.5mm	STD	Plano
8.20mm	14.5mm	FLT2	Plano
8.60mm	14.5mm	STP2	Plano

Risk-Free Warranty

KeraSoft® Thin quarterly replacement silicone hydrogel lenses are available in single lens vials with an optional warranty that provides 2 exchanges per eye and full cancellation privileges (less shipping & handling) within 90 days from initial order. After successfully fitting a patient using the warranty, save by ordering non-warranted lenses in the quantity appropriate to the prescribed replacement schedule.

Shipping, Credit & Return Policies

Standard production time for KeraSoft® Thin lenses is 4 business days. Lenses will be shipped via Ground service unless other traceable option is requested. Warranted KeraSoft® Thin lenses being returned for exchange or credit must meet the following guidelines:

- returned in original glass vial with original label
- returned within specified warranty period
- accompanied by a completed Request for Credit Form (provided with each order) exact reason for return must be specified and original invoice number & order date must be included

Non-warranted lenses are non-returnable/non-refundable. Shipping fees and proof of delivery on all returned product is the responsibility of the customer. Credit is not issued for shipping & handling charges.

