

Advanced Designs. Simple Success.

MetroSoft Soft Lenses

Solid performance and reliability are trademarks of the MetroSoft line of lenses. They are now available in Polymacon 38% water content, Hioxifilcon B 49% water content, Definitive™ Silicone Hydrogel material, 60Dk, 74% water content these materials provide durability and high visual acuity, optical clarity, stability and all day comfort. These lenses are custom made on CNC lathing systems that provide the consistency and accuracy necessary to meet the needs of your soft contact lens patients. From the highly astigmatic or the astigmatic presbyope, the MetroSoft line of lenses can provide patient retention though the consistency and quality of our manufacturing process.

METROSOFT TORIC

Base Curve: 8.4, 8.7, 9.0 (polymacon/definitive)
 8.1, 8.4, 8.7 (hioxifilcon B)
 Power: -12.00 to +12.00
 Cylinder: -.75 to -8.00 in .25 steps
 Axis: 5 to 180 in 5 degree steps
 Diameter: 14.0 (polymacon)
 14.2 (definitive/hioxifilcon B)
 Prism: one diopter
 Color: Clear or Light Blue

MetroFocal MULTIFOCAL SOFT LENS

Base Curve: 8.3, 8.6, 8.9 (polymacon/definitive)
 8.1, 8.4, 8.7 (hioxifilcon B)
 Power: -12.00 to +12.00
 Diameter: 14.0 (polymacon)
 14.2 (definitive/hioxifilcon B)
 Add Power: up to +2.50
 Color: Clear or Light Blue

MetroFocal Toric MULTIFOCAL SOFT LENS

Base Curve: 8.3, 8.6, 8.9 (polymacon/definitive)
 8.1, 8.4, 8.7 (hioxifilcon B)
 Power: -12.00 to +12.00
 Cylinder: -.75 to -8.00
 Axis: 5 to 180 in 5 degree steps
 Diameter: 14.0 (polymacon)
 14.2 (definitive/hioxifilcon B)
 Add Power: up to +2.50
 Prism: one diopter
 Color: Clear or Light Blue

MetroSoft II

Base Curve: 8.3, 8.6, 8.9 (polymacon/definitive)
 8.1, 8.4, 8.7 (hioxifilcon B)
 Power: -12.00 to +12.00
 Diameter: 13.5 (polymacon)
 14.2 (definitive/hioxifilcon B)
 Color: Clear or Light Blue

MetroTint Available only in polymacon

Base Curve: 8.3mm, 8.6mm, 8.9mm
 Power: -12.00 to +12.00
 Diameter: 13.5mm
 Color: Aqua, Blue, Green

Extended parameters available